

## Triumph IP12 Infusion Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 14 preinstalled IV lines to include Baxter, Covetrus A Series, Covetrus B Series, EcoMed B Series, ICU Medical, Mvet A Series, Mvet B Series, Pivetal A Series, PivetalB Series, VetOne A Series, VetOne B Series, Zoetis A Series, and Zoetis B Series.

### Product Downloads:

- IP12 Infusion Pump Spec Sheet
- IP12 Infusion Pump Quick Start Guide
- IP12 Infusion Pump User Manual

### Basic Infusion Startup Process - Triumph IP12 Infusion Pump

### Infusion Modes Training Video for Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump



medical services, llc

## **How To Customize Settings for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump**

## **Sequence Mode Programming for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump**

## **Common User Errors for the Triumph IP12 Infusion Pump**