

Triumph SP12 Syringe Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 8 preinstalled Syringes to include Monoject, BD, Terumo, VetOne, Covetrus, Pivotal, EcoMed, and Exel

Product Downloads:

- SP12 Syringe Pump Spec Sheet
- SP12 Syringe Pump Quick Start Guide
- SP12 Syringe Pump User Manual

Basic Infusion Startup Process - Triumph SP12 Syringe Pump

Infusion Modes Training Video for Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump

How To Customize Settings for the Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump



How to Correctly Load a Syringe in the Triumph SP12 Syringe Pump

Sequence Mode Programming for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump