

## Triumph SP12 Syringe Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 8 preinstalled Syringes to include Monoject, BD, Terumo, VetOne, Covetrus, Pivotal, EcoMed, and Exel

### Product Downloads:

- SP12 Syringe Pump Spec Sheet
- SP12 Syringe Pump Quick Start Guide
- SP12 Syringe Pump User Manual
- Acceptable and Current Syringe Brands and Sizes For The Triumph SP12\*

\*Using a syringe brand and/or size not listed or currently calibrated in the Triumph SP12 can void your warranty.

### Basic Infusion Startup Process - Triumph SP12 Syringe Pump

### Infusion Modes Training Video for Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump



## **How To Customize Settings for the Triumph SP12 Syringe Pump and Triumph IP12 Infusion Pump**

### **How to Correctly Load a Syringe in the Triumph SP12 Syringe Pump**

### **Sequence Mode Programming for the Triumph IP12 Infusion Pump and Triumph SP12 Syringe Pump**