

Triumph TRI-TFW26 Vet Fluid Warmer

The Triumph TRI-TFW26 Vet Fluid Warmer is a 'dry heater', which does not directly contact the heated liquid. The heat is transferred from the warming plates through the infusion line. The user can adjust the set temperature in the range of 30°C - 50°C according to actual needs. The warmer maintains a constant temperature state after reaching the target temperature.

- **Compact size, better appearance**
- **More accurate temperature control**
- **Easy to use and operate**
- **30° C - 50° C temperature range, adjustable**
- **Two mounting methods: IV stand or with the provided lanyard and bracket**
- **Over-temperature protection, five alarms including high and low temperature alarm, safer and more reliable.**

Product Downloads:

- [Triumph Vet Fluid Warmer Spec Sheet](#)
- [Triumph Vet Fluid Warmer User Manual](#)