

Triumph VIP Infusion Pump

- Large 4.3 inch water resistant color touch screen
- Dual Battery for providing over 9 hours of battery life
- Enhanced safety features including full screen color coded alarms
- Rate Mode, Time Mode, Body-Weight Mode, Drug Library Mode, and Drip Mode
- Fully customizable drug library with password protection
- Folding carry handle providing easier transport and space saving stackable design
- Universal pole clamp enabling cage mounting
- 14 preinstalled IV lines to include Baxter, Covetrus A Series, Covetrus B Series, EcoMed B Series, ICU Medical, Mvet A Series, Mvet B Series, Pivotal A Series, Pivotal B Series, VetOne A Series, VetOne B Series, Zoetis A Series, and Zoetis B Series.

Product Downloads:

- VIP Infusion Pump Spec Sheet
- VIP Infusion Pump Quick Start Guide
- VIP Infusion Pump User Manual

Triumph VIP Infusion Pump Basic Setup

Triumph VIP Infusion Pump Features & Settings

Triumph VIP Infusion Pump Drug Library



Triumph VIP Infusion Pump Features During Infusion

Triumph VIP Infusion Pump Extended Warranty Offer

Triumph VIP Pump Common User Error with Power Cord