

# Vetland Landmark VSA-2100

The small, stable, and versatile anesthesia solution. This device meets all of your anesthesia needs with a small footprint. It includes a 5-leg base with smooth casters and is very lightweight. These features promise one of the most stable products available for any use.

## Features

- Neat, balanced, compact design
- Ships with a Drager Isoflurane Funnel Fill Vaporizer
- Pop-off valve with thumb occlusion
- Proven human quality FDA type valves
- Best heat and moisture conserving coaxial absorber
- 1030cc absorber is better
  - Less trapped volume
  - Less chance of “channeling”
  - Quicker anesthetic adjustments
- Quick disconnect fresh gas outlet for non-rebreathing applications
- Mayo tray standard
- Back cover protects all tubing connections
- 0-4 lpm oxygen flow control with protective stop
- Corrosion proof construction
- 5-leg base with lifetime caster replacements

## **Starter Kit (Required with new rental)**

- 1 Adult 22mm F circuit with 2L Breathing Bag
- 1 Pediatric 15mm wye circuit with 1 L Breathing Bag
- 1 Bag of CO2 absorbent (Litholyme)
- 1 6-foot, 19mm, scavenger tubing to connect to F/Air
- 1 19mm to 22mm adaptor for F/Air to mate up to Scavenger tubing

\* Green listed items are consumable accessories, additional accessories are purchased through your MWI rep.