

# **Zevex EnteraLite Infinity**

Dimensions: 1.95"H x 5.65"W x 4.05"D

Weight: 14.4 oz

Flow Rate Range:

0.1-9.9 ml/hr in 0.1 ml/hr increments

0-600 ml/hr in 1 ml/hr increments

Dose Range:

0.1-9.9 ml in 0.1 ml increments

10-999 ml in 1 ml increments

1000-3000 ml in 10 ml increments

Accuracy:  $\pm 5\%$

Battery: Rechargeable Lithium Ion

Battery Life: 24 hours at 125 ml/hr

Recharge Time: Approximately 6 hours